



BIOFECT

RENGÖRINGSMEDEL

BIOFECT DETERGENT

MATERIAL CARE

Disinfectant
detergent

DISINFECTANT

APPLICATION

To be used for sanitising various environments to inhibit the spread of infection. Also functions as a detergent.

Eclipse Biofarmab Biofect is a highly concentrated, alkaline liquid cleaner that offers both cleaning and disinfecting properties. The product is effective at low temperatures, is without risk to animals and humans and has no negative effects on metals.

Suitable for a wide range of applications, Biofect is ideal for cleaning animal stables, boxes, horse trailers, veterinary clinics, and other animal transport vehicles, as well as any space requiring stringent infection control measures. Additionally, it can be used as a detergent in washing machines for laundering blankets, clothing, and various types of equipment.

The active substance falls under the category of quaternary ammonium compounds. It has been proven effective against a variety of pathogens, including:

Gram-positive bacteria (duration of action 15 minutes)

- *Staphylococcus Aureus*
- *Listeria monocytogenes*

Gram-negative bacteria (duration of action 5 minutes)

- EHEC (*Enterohaemorrhagic Escherichia coli*)
- *Salmonella Typhimurium*
- *Campylobacter jejuni*

Additionally, it is effective against:

- Virus
- Fungi/yeast *Candida Albicans*
- Mold/mold spores *Aspergillus Niger*

INSTRUCTIONS FOR USE

For surface cleaning:

1. Remove loose dirt.
2. Rinse all surfaces with water; scrub if needed.
3. For regular cleaning: Mix 1 dl with 10 l of water.
- For tough cleaning: Mix 2.5 dl with 10 l of water.
4. Apply the solution to all surfaces and let it sit for approximately 20 minutes.
5. Rinse thoroughly with clean water and allow the surface to dry.

For washing machine:

Heavily soiled laundry: Use 100 ml.
Normally soiled laundry: Use 50 ml.

Use protective safety glasses and gloves when handling concentrated product to avoid contact with eyes and skin.

SIZES AVAILABLE

0550 2.5 l

COMPOSITION

Alkylbenzyltrimethylammonium chloride, 5%—15%
Cationic Surfactants.

BIOCIDALLY ACTIVE SUBSTANCE

7.5g Alkyl (C12-16) Dimethylbenzylammonium Chloride (ADBAC/BKC (C12C16)) per 100 g Liquid Concentrate Disinfectant.

SHELF LIFE

24 months.

